i-BONDING® LC^N

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EN

INDICATIONS

- **Used to bond dentin and enamel to polymer-based restorative materials, compomers, ceramic surfaces to improve retention of a restoration.**

DESCRIPTION

i-BONDING LC^N is a 5th generation light curing single component one-step adhesive based on acetone and resins HEMA and Bis-GMA. Is suitable for use with wet bonding technique. Etching of the surfaces is required.

INSTRUCTION FOR USE

NOTE! Do not use if primary package is broken. It is recommended to use cofferdam during application of the product.

1. Prepare cavity as always. Cleaning with oil-free prophylaxis paste, such as i-FASTE is recommended. For deep cavities use calcium hydroxide liner or glass ionomer based cement, such as i-BAS.
2. Apply layer of etch, such as i-GEL^N on all surface to be etched. Leave etch in place for 15 seconds (dentine), 30 seconds (enamel). Rinse with water and dry with air. Avoid over drying dentin. After drying, the surface must be chalky and do not contaminated prior to applying the bonding resin. If the surface becomes contaminated with saliva, rinse and dry again and re-etch if necessary.
3. Shake bottle before each use.
4. Delicate push on bottle and dispense small amount of adhesive into cup. Do not force bottle.
5. Using a fully saturated micro applicator apply 2-3 layers of adhesive onto etched surface.
6. Dry with air for 10-20 seconds to dry all solvent from surface. If surface is not glossy, apply additional layers of adhesive and dry.
7. Curing with light for 20 seconds. Use any polymerization lamp with light intensity $\geq 700\text{mW/cm}^2$ and the wavelength range of 400-500nm (blue light).
8. Apply the layer of light curing restorative material, such as i-LIGHT^N, i-XCITE LC^N or i-FLOW^N into cavity and finish restoration.

NOTE! After the desired amount of adhesive extruded, immediately close the bottle.

WARNINGS AND PRECAUTIONS!

NOTE! i-BONDING LC^N is highly flammable liquid. Do not use i-BONDING LC^N if patient is allergic to product or to any of the ingredients. Not to be used on exposed pulp.

Hazard and precautionary statements:

Uncured product causes skin, serious eye irritation, may cause respiratory irritation. IF ON SKIN: Wash with plenty of soap and water. If skin irritation occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a POISON CENTRE or doctor/physician if you feel unwell.

Prevention:

Avoid breathing vapours. Wash hands thoroughly after handling. Use only in a well-ventilated area. Wear protective gloves/protective clothing/eye protection/face protection.

Interaction with other substances:

i-BONDING LC^N should not be used with products containing eugenol.

STORAGE

Keep container tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight and heat sources. Do not freeze the product. The batch number should be quoted in all correspondence. See packaging for batch and expiry date. Do not use after the expiry date. Dispose of contents/container to as required by national regulatory requirements. Keep out of the reach of children!

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

PACKAGING

REF IBLTP 5ml bottle, 5pcs micro applicators, cup

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