

Disposable Micro Applicators

CE

EN Disposable micro applicators



INSTRUCTIONS FOR USE

EN

For application of sealants, adhesives, etchants, cavity liners and other different solutions. Bends into any angle. Three head sizes: regular (2.5mm), fine (2.00mm), ultrafine (1.5mm).

Choose size of the disposable micro applicators head depending on the substance and area accessibility. Immerse disposable micro applicators head into the substance and then apply onto the selected surface.

WARNINGS AND PRECAUTIONS!

For single use only.

STORAGE

Keep tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight, sources of ignition. Do not use after the expiry date. Dispose of contents/container to as required by national regulatory requirements. Keep out of the reach of children!

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

PACKAGING

REF **IB05R** Regular size, colours: blue and light blue, 100pcs.

REF **IB05F** Fine size, colours: green and light green, 100pcs.

REF **IB05U** Ultrafine size, colours: lavender and light lavender, 100pcs.

MANUFACTURER:



i-dental

UAB „Medicinos linija“
Karaliaučiaus str. 29, LT-78374 Šiauliai, Lithuania
Tel. +370 41 553 553, fax +370 41 553 551
dental@i-dental.lt, www.i-dental.lt

Last revised: 2014-11