

Mixing pads



EN Dental mixing pads

LT Padėkliukai odontologinių medžiagų ir preparatų maišymui



INSTRUCTIONS FOR USE

EN

Tear-off, disposable coated sheets for the easy mixing of materials by preventing absorption into the surface. Place the materials to be mixed on the mixing pad. Mix with suitable instrument. After mixed material is used tear-off and dispose the sheet.

WARNINGS AND PRECAUTIONS!

For single use only.

STORAGE

Keep tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight, sources of ignition. Do not use after the expiry date. Dispose of contents/container to as required by national regulatory requirements. Keep out of the reach of children!

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

PACKAGING

REF I B04 5.1x5.1cm, 20 sheets

REF I B041 7.6x7.6cm, 50 sheets

NAUDOJIMO INSTRUKCIJA

LT

Nuplėšiami, vienkartiniai lapeliai, skirti medžiagų maišymui. Paviršius padengtas, kad medžiagos neįsigertų. Maišymui skirtą medžiagą padėkite ant padėkliuko. Maišykite su tam skirtu instrumentu. Sumaišę medžiagas panaudotą lapelį nuplėškite ir išmeskite.

ĮSPĖJIMAI IR ATSARGUMO PRIEMONĖS!

Tik vienkartiniam naudojimui.

LAIKYMAS

Laikykite sandariai uždaryta, sausoje, gerai vėdinamoje vietoje, 4-28°C temperatūroje. Saugokite nuo tiesioginių saulės spindulių ir šilumos šaltinių. Nenaudokite produkto pasibaigus jo galiojimo terminui. Turinį ar talpyklą išmeskite laikydamiesi nacionalinių norminių reikalavimų. Laikykite vaikams nepasiekiamoje vietoje!

GAMINTOJO ATSAKOMYBĖ

Mūsų produktas sukurtas profesionaliam naudojimui odontologijoje. Mes negalime kontroliuoti kaip šis produktas naudojamas, todėl už tinkamą jo panaudojimą atsako vartotojas. Žinoma, mes užtikriname, produkto kokybės atitikimą taikomiems standartams.

PAKUOTĖ

REF I B04 5.1x5.1cm, 20 lapelių

REF I B041 7.6x7.6cm, 50 lapelių

MANUFACTURER/ GAMINTOJAS:



i-dental

UAB „Medicinos linija“
Karaliaučiaus str. 29,
LT-78348 Šiauliai, Lithuania
Tel. +370 41 553 553
info@i-dental.lt, www.i-dental.lt

Last revised/ Paskutinė peržiūra: 2018-02